

Diethylstilbestrol (cis/trans mix) (ring-3,3',5,5'-diethyl-1,1,1',1'-D8,98%) 100 µg/mL in Dioxane

Art. ID CIL-DLM-170-D-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1165 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Diethylstilbestrol (cis /trans mix) (ring-3,3', 5,5'-diethyl-1,1,1',1'-D8,98%)	[91318-10-4]	100	µg/mL		