

2,4-Dinitrophenol (ring-D3,98%) (contains 0.35 mg/mL deuterium) 1000 µg/mL in Deuterium oxide

Art. ID CIL-DLM-299-10
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4-Dinitrophenol (ring -D3,98%)	[93951-77-0]	1000	µg/mL		