

**3,3',4,4'-Tetrachlorobiphenyl (IUPAC #77) (D6,98%) 40 µg/mL in Nonane**

Art. ID CIL-DLM-3063-3  
Unit 3 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 97%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,4'-Tetrachlorobiphenyl (IUPAC #77) (D6, 98%)	[93952-23-9]	40	µg/mL		