

**Dibenzo(a,i)pyrene (D14,98%) 200 µg/mL in Toluene-d8**

Art. ID CIL-DLM-3740-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dibenzo(a,i)pyrene (D14,98%)	[158776-07-9]	200	µg/mL		