

**p-Terphenyl (D14,98%) 200 µg/mL in Isooctane**

Art. ID                    CIL-DLM-382-1.2  
Unit                        1,2 mL  
Deliverydetails        ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

**Description**

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	p-Terphenyl D14	[1718-51-0]	200	µg/mL		