

Nitroacenphthene (D9,98 %) 50 µg/mL in Toluene-d8

Art. ID CIL-DLM-3836-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Nitroacenphthene (D9,98 %)	[2483831-38-3]	50	µg/mL		