

### **Fluorene, 2-nitro (D9,98%) 50 µg/mL in Toluene-D8**

Art. ID CIL-DLM-3837-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

#### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Fluorene, 2-nitro (D9,9 8%)	[128008-87-7]	50	µg/mL		