

6-Nitrochrysen (D11,98 %) 50 µg/mL in Toluene-D8

Art. ID CIL-DLM-3839-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Description

Chemical purity: 95%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	6-Nitrochrysen (D11,98 %)	[203805-92-9]	50	µg/mL		