

Dibenzo(a,j)acridine (D13,98%) 50 µg/mL in Toluene-d8

Art. ID CIL-DLM-3843-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dibenzo(a,j)acridine (D 13,98%)	[2483831-41-8]	50	µg/mL		