

Genistein (3',5',6,8-D4,95%) 100 ± 10 µg/mL in Acetonitrile

Art. ID CIL-DLM-4460-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Genistein (3',5',6,8-D4 ,95%)	[136466-47-2]	100	µg/mL		