

### **Phosmet (dimethyl-D6,98%) 100 µg/mL in Acetonitrile**

Art. ID CIL-DLM-4667-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

#### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosmet (dimethyl-D6,98 %)	[2083623-41-8]	100	µg/mL		