

### **3-Nitrofluoranthene (D9,98%) 50 µg/mL in Toluene-D8**

Art. ID CIL-DLM-4711-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

#### Description

Chemical purity: 87% (approx. 13% OTHER ISOMER)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-Nitrofluoranthene (D9,98%)	[350820-11-0]	50	µg/mL		