

## **2,4,6-Trichloroanisole (D5,98%) 1 mg/mL in Methanol-D4**

Art. ID CIL-DLM-6083-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4,6-Trichloroanisole (D5,98%)	[352439-08-8]	1	mg/mL		