

Tris (2-chloroisopropyl) phosphate (D18,98%) 100 µg/mL in Acetonitrile

Art. ID CIL-DLM-9317-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tris (2-chloroisopropyl) phosphate (D18,98%)	[13674-84-5]/un labelled	100	µg/mL		