

Coplanar PCB Mixture

Art. ID CIL-EC-4070
Unit 3 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 77	3,3',4,4'-Tetrachlorobiphenyl (13C12,99%) (IUPAC #77)	[105600-23-5]	5	µg/mL		
IUPAC# 126	3,3',4,4',5-Pentachlorobiphenyl (13C12,99%) (IUPAC #126)	[208263-65-4]	5	µg/mL		
IUPAC# 169	3,3',4,4',5,5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #169)	[208263-70-1]	5	µg/mL		