

## Coplanar PCB Standard Solution

Art. ID CIL-EC-4187  
Unit 3 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 77	3,3',4,4'-Tetrachlorobiphenyl (13C12,99%) (IUPAC #77)	[105600-23-5]	1	µg/mL		
IUPAC# 81	3,4,4',5-Tetrachlorobiphenyl (13C12,99%) (IUPAC #81)	[208461-24-9]	1	µg/mL		
IUPAC# 126	3,3',4,4',5-Pentachlorobiphenyl (13C12,99%) (IUPAC #126)	[208263-65-4]	1	µg/mL		
IUPAC# 169	3,3',4,4',5,5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #169)	[208263-70-1]	1	µg/mL		