

## **Method 1668A/B Labelled Injection Internal Standard Solution**

Art. ID CIL-EC-4979  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 9L	2,5-Dichlorobiphenyl (13C12,99%) (IUPAC #9)	[250694-89-4]	5	µg/mL		
IUPAC# 52L	2,2',5,5'-Tetrachlorobiphenyl (13C12,99%) (IUPAC #52)	[208263-80-3]	5	µg/mL		
IUPAC# 101L	2,2',4,5,5'-Pentachlorobiphenyl (13C12,99%) (IUPAC #101)	[104130-39-4]	5	µg/mL		
IUPAC# 138L	2,2',3,4,4',5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #138)	[208263-66-5]	5	µg/mL		
IUPAC# 194L	2,2',3,3',4,4',5,5'-Octachlorobiphenyl (13C12,99%) (IUPAC #194)	[208263-74-5]	5	µg/mL		