

## **PCB Mixture (PCB-70/111/138/170)**

Art. ID CIL-EC-5163  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 70	2,3',4',5-Tetrachlorobiphenyl (13C12,99%) (IUPAC #70)	[208263-81-4]	1000	ng/mL		
IUPAC# 111	2,3,3',5,5'-Pentachlorobiphenyl (13C12,99%) (IUPAC #111)	[235416-29-2]	1000	ng/mL		
IUPAC# 138	2,2',3,4,4',5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #138)	[208263-66-5]	1000	ng/mL		
IUPAC# 170	2,2',3,3',4,4',5'-Heptachlorobiphenyl (13C12,99%) (IUPAC #170)	[35065-30-6]/unlabelled	1000	ng/mL		