

Modified JiS PCB Sampling Spike

Art. ID CIL-EC-5326
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 79	3,3',4,5'-Tetrachlorobi phenyl (13C12,99%) (IUP AC #79)	[41464-48-6]/un labelled	50	ng/mL		