

## **EN-1948-4 PCB sampling standard**

Art. ID	CIL-EC-5370
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 60	2,3,4,4'-Tetrachlorobiphenyl (13C12,99%) (IUPAC #60)	[n/a]	100	ng/mL		
IUPAC# 127	3,3',4,5,5'-Pentachlorobiphenyl (13C12,99%) (IUPAC #127)	[n/a]	100	ng/mL		
IUPAC# 159	2,3,3',4,5,5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #159)	[n/a]	100	ng/mL		