

**EN-1948-4 PCB sampling standard (1/10 concentration)**

Art. ID CIL-EC-5370-1/10X-10  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 60	2,3,4,4'-Tetrachlorobiphenyl (13C12,99%) (IUPAC #60)	[n/a]	10	ng/mL		
IUPAC# 127	3,3',4,5,5'-Pentachlorobiphenyl (13C12,99%) (IUPAC #127)	[n/a]	10	ng/mL		
IUPAC# 159	2,3,3',4,5,5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #159)	[n/a]	10	ng/mL		