

## **EN-1948-4 PCB recovery standard**

Art. ID	CIL-EC-5371
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 70	2,3',4',5-Tetrachlorobiphenyl (13C12,99%) (IUPAC #70)	[208263-81-4]	100	ng/mL		
IUPAC# 111	2,3,3',5,5'-Pentachlorobiphenyl (13C12,99%) (IUPAC #111)	[235416-29-2]	100	ng/mL		
IUPAC# 170	2,2',3,3',4,4',5-Heptachlorobiphenyl (13C12,99%) (IUPAC #170)	[35065-30-6]/unlabelled	100	ng/mL		