

## **EN-1948-4 Marker PCB extraction standard**

Art. ID	CIL-EC-5379
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 28	2,4,4'-Trichlorobiphenyl (13C12,99%) (IUPAC #28)	[208263-76-7]	1000	ng/mL		
IUPAC# 52	2,2',5,5'-Tetrachlorobiphenyl (13C12,99%) (IUPAC #52)	[208263-80-3]	1000	ng/mL		
IUPAC# 101	2,2',4,5,5'-Pentachlorobiphenyl (13C12,99%) (IUPAC #101)	[104130-39-4]	1000	ng/mL		
IUPAC# 138	2,2',3,4,4',5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #138)	[208263-66-5]	1000	ng/mL		
IUPAC# 153	2,2',4,4',5,5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #153)	[185376-58-3]	1000	ng/mL		
IUPAC# 180	2,2',3,4,4',5,5'-Heptachlorobiphenyl (13C12,99%) (IUPAC #180)	[208263-72-3]	1000	ng/mL		