

## **EN-1948-4 Marker PCB extraction standard**

Art. ID	CIL-EC-5379-1/10X-10
Unit	10 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 28	2,4,4'-Trichlorobipheny l (13C12,99%) (IUPAC #2 8)	[208263-76-7]	100	ng/mL		
IUPAC# 52	2,2',5,5'-Tetrachlorobi phenyl (13C12,99%) (IUP AC #52)	[208263-80-3]	100	ng/mL		
IUPAC# 101	2,2',4,5,5'-Pentachloro biphenyl (13C12,99%) (I UPAC #101)	[104130-39-4]	100	ng/mL		
IUPAC# 138	2,2',3,4,4',5'-Hexachlo robiphenyl (13C12,99%) (IUPAC #138)	[208263-66-5]	100	ng/mL		
IUPAC# 153	2,2',4,4',5,5'-Hexachlo robiphenyl (13C12,99%) (IUPAC #153)	[185376-58-3]	100	ng/mL		
IUPAC# 180	2,2',3,4,4',5,5'-Heptac hlorobiphenyl (13C12,99 %) (IUPAC #180)	[208263-72-3]	100	ng/mL		