

2,3,3',5-Tetrachlorobiphenyl 35 µg/mL in Isooctane

Art. ID CIL-PCB-58
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------------|--------------|---------------------|-------|--------|--------|
| | 2,3,3',5-Tetrachlorobiphenyl (PCB 57) | [70424-67-8] | 35 | µg/mL | | |