

## **2,3,3',6-Tetrachlorobiphenyl 35 µg/mL in Isooctane**

Art. ID CIL-PCB-59  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

### Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter                     | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|---------------------------------------|--------------|---------------------|-------|--------|--------|
|                  | 2,3,3',6-Tetrachlorobiphenyl (PCB 59) | [74472-33-6] | 35                  | µg/mL |        |        |