

## **2,3,4,4'-Tetrachlorobiphenyl 100 µg/mL in Isooctane**

|                 |  |
|-----------------|--|
| Art. ID         | CIL-PCB-60-CS  |
| Unit            | 1,2 mL   |
| Deliverydetails | ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2 |

### Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter                     | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|---------------------------------------|--------------|---------------------|-------|--------|--------|
|                  | 2,3,4,4'-Tetrachlorobiphenyl (PCB 60) | [33025-41-1] | 100                 | µg/mL |        |        |