

3,3',4,5-Tetrachlorobiphenyl 100 µg/mL in Isooctane

Art. ID	CIL-PCB-78-CS
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,5-Tetrachlorobiphenyl (PCB 78)	[70362-49-1]	100	µg/mL		