

### **3,3',4,5'-Tetrachlorobiphenyl 100 µg/mL in Isooctane**

Art. ID CIL-PCB-79-CS  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

#### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,5'-Tetrachlorobiphenyl (PCB 79)	[41464-48-6]	100	µg/mL		