

3,4,4',5-Tetrachlorobiphenyl 100 µg/mL in Isooctane

Art. ID	CIL-PCB-81-CS
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,4,4',5-Tetrachlorobiphenyl (PCB 81)	[70362-50-4]	100	µg/mL		