

2,2',3,4,4'-Pentachlorobiphenyl 100 µg/mL in Isooctane

| | |
|-----------------|--|
| Art. ID | CIL-PCB-85-CS |
| Unit | 1,2 mL |
| Deliverydetails | ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2 |

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|--------------|---------------------|-------|--------|--------|
| | 2,2',3,4,4'-Pentachloro biphenyl (PCB 85) | [65510-45-4] | 100 | µg/mL | | |