

2,2',3,4,4'-Pentachlorobiphenyl 100 µg/mL in Isooctane

Art. ID	CIL-PCB-85-CS
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,4,4'-Pentachloro biphenyl (PCB 85)	[65510-45-4]	100	µg/mL		