

2,2',3,4,5'-Pentachlorobiphenyl 100 µg/mL in Isooctane

Art. ID	CIL-PCB-87-CS
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,4,5'-Pentachloro biphenyl (PCB 87)	[38380-02-8]	100	µg/mL		