

Modified JIS PCB alternate B syringe spike

Art. ID CIL-EC-5420
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|-------------------------|---------------------|-------|--------|--------|
| IUPAC# 70L | 2,3',4',5-Tetrachlorobiphenyl (13C12,99%) (IUPAC #70) | [208263-81-4] | 50 | ng/mL | | |
| IUPAC# 111L | 2,3,3',5,5'-Pentachlorobiphenyl (13C12,99%) (IUPAC #111) | [235416-29-2] | 50 | ng/mL | | |
| IUPAC# 138L | 2,2',3,4,4',5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #138) | [208263-66-5] | 50 | ng/mL | | |
| IUPAC# 170L | 2,2',3,3',4,4',5'-Heptachlorobiphenyl (13C12,99%) (IUPAC #170) | [35065-30-6]/unlabelled | 50 | ng/mL | | |