

## **DL-PCB RH12 Syringe Spike**

Art. ID CIL-EC-5423  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 70L	2,3',4',5-Tetrachlorobiphenyl (13C12,99%) (IUPAC #70)	[208263-81-4]	100	ng/mL		
IUPAC# 127L	3,3',4,5,5'-Pentachlorobiphenyl (13C12,99%) (IUPAC #127)	[n/a]	100	ng/mL		
IUPAC# 128L	2,2',3,3',4,4'-Hexachlorobiphenyl (13C12,99%) (IUPAC #128)	[38380-07-3]/unlabelled	100	ng/mL		
IUPAC# 170L	2,2',3,3',4,4',5-Heptachlorobiphenyl (13C12,99%) (IUPAC #170)	[35065-30-6]/unlabelled	100	ng/mL		