

## **Rapid PCB Screening Calibration Solutions CS4 (unlabelled/<sup>13</sup>C12,99%)**

Art. ID CIL-EC-5448-CS4  
Unit 0,2 mL  
Deliverydetails De Minimis UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
CS4, IUPAC# 28	2,4,4'-Trichlorobipheny I (PCB 28)	[7012-37-5]	50	ng/mL		
CS4, IUPAC# 18	2,2',5-Trichlorobipheny I (PCB 18)	[37680-65-2]	50	ng/mL		
CS4, IUPAC# 44	2,2',3,5'-Tetrachlorobi phenyl (PCB 44)	[41464-39-5]	50	ng/mL		
CS4, IUPAC# 70	2,3',4',5-Tetrachlorobi phenyl (PCB 70)	[32598-11-1]	50	ng/mL		
CS4, IUPAC# 52	2,2',5,5'-Tetrachlorobi phenyl (PCB 52)	[35693-99-3]	50	ng/mL		
CS4, IUPAC# 101	2,2',4,5,5'-Pentachloro biphenyl (PCB 101)	[37680-73-2]	50	ng/mL		
CS4, IUPAC# 110	2,3,3',4',6-Pentachloro biphenyl (PCB 110)	[38380-03-9]	50	ng/mL		
CS4, IUPAC# 118	2,3',4,4',5-Pentachloro biphenyl (PCB 118)	[31508-00-6]	50	ng/mL		
CS4, IUPAC# 149	2,2',3,4',5',6-Hexachlo robiphenyl (PCB 149)	[38380-04-0]	50	ng/mL		
CS4, IUPAC# 138	2,2',3,4,4',5'-Hexachlo robiphenyl (PCB 138)	[35065-28-2]	50	ng/mL		
CS4, IUPAC# 153	2,2',4,4',5,5'-Hexachlo robiphenyl (PCB 153)	[35065-27-1]	50	ng/mL		
CS4, IUPAC# 180	2,2',3,4,4',5,5'-Heptac hlorobiphenyl (PCB 180)	[35065-29-3]	50	ng/mL		
CS4, IUPAC# 187	2,2',3,4',5,5',6-Heptac hlorobiphenyl (PCB 187)	[52663-68-0]	50	ng/mL		
CS4, IUPAC# 28L	2,4,4'-Trichlorobipheny I ( <sup>13</sup> C12,99%) (IUPAC #2 8)	[208263-76-7]	10	ng/mL		

CS4, IUPAC# 70L	2,3',4',5'-Tetrachlorobiphenyl (13C12,99%) (IUPAC #70)	[208263-81-4]	10	ng/mL
CS4, IUPAC# 52L	2,2',5,5'-Tetrachlorobiphenyl (13C12,99%) (IUPAC #52)	[208263-80-3]	10	ng/mL
CS4, IUPAC# 101L	2,2',4,5,5'-Pentachlorobiphenyl (13C12,99%) (IUPAC #101)	[104130-39-4]	10	ng/mL
CS4, IUPAC# 118L	2,3',4,4',5'-Pentachlorobiphenyl (13C12,99%) (IUPAC #118)	[104130-40-7]	10	ng/mL
CS4, IUPAC# 138L	2,2',3,4,4',5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #138)	[208263-66-5]	10	ng/mL
CS4, IUPAC# 153L	2,2',4,4',5,5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #153)	[185376-58-3]	10	ng/mL
CS4, IUPAC# 141L	2,2',3,4,5,5'-Hexachlorobiphenyl (13C12,99%) (IUPAC #141)	[n/a]	10	ng/mL
CS4, IUPAC# 180L	2,2',3,4,4',5,5'-Heptachlorobiphenyl (13C12,99%) (IUPAC #180)	[208263-72-3]	10	ng/mL