

Rapid PCB Screening Syringe Spike (13C12,99%)

Art. ID	CIL-EC-5450
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
IUPAC# 70L	2,3',4',5-Tetrachlorobi phenyl (13C12,99%) (IUP AC #70)	[208263-81-4]	2000	ng/mL		
IUPAC# 70L	2,3',4,4',5-Pentachloro biphenyl (13C12,99%) (I UPAC #118)	[104130-40-7]	2000	ng/mL		
IUPAC# 70L	2,2',3,4,5,5'-Hexachlor obiphenyl (13C12,99%) (I UPAC #141)	[n/a]	2000	ng/mL		