

## **Chlorinated Internal Standard Stock Solution 1 (unlabelled)**

Art. ID CIL-EC-5477  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tetrachloro-m-xylene	[877-09-8]	25	µg/mL		
IUPAC# 209	2,2',3,3',4,4',5,5',6,6'-Decachlorobiphenyl (P CB 209)	[2051-24-3]	25	µg/mL		