

PXB Calibration Solution [CS3]

Art. ID	CIL-ECB-5390-CS3
Unit	0,2 mL
Deliverydetails	De Minimis UN1920 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
CS3	4'-Bromo-3,3',4,5-tetra chlorobiphenyl	[n/a]	50	ng/mL		
CS3	4'-Bromo-2,3',4,5-tetra chlorobiphenyl	[n/a]	50	ng/mL		
CS3	4'-Bromo-2,3,3',4-Tetra chlorobiphenyl		50	ng/mL		
CS3	4'-bromo-2,3,3',4,5-pen tachlorobiphenyl	[n/a]	50	ng/mL		
CS3	4'-Bromo-3,3',4,5,5'-pe ntachlorobiphenyl	[n/a]	50	ng/mL		
CS3	3,4-Dichloro-3',4',5'-t ribromobiphenyl	[n/a]	100	ng/mL		
CS3	4'-Bromo-3,3',4,5-tetra chlorobiphenyl (13C12,9 9%)		100	ng/mL		
CS3	4'-Bromo-2,3',4,5-tetra chlorobiphenyl (13C12,9 9%)		100	ng/mL		
CS3	4'-Bromo-2,3,3',4-tetra chlorobiphenyl (13C12,9 9%)		100	ng/mL		
CS3	4'-Bromo-2,3,3',4,5-pen tachlorobiphenyl (13C12 ,99%)		100	ng/mL		
CS3	4'-Bromo-3,3',4,5,5'-pe ntachlorobiphenyl (13C1 2,99%)		100	ng/mL		
CS3	3,4-Dichloro-3',4',5'-t ribromobiphenyl (13C12,		200	ng/mL		

99%)			
CS3	2,2',3,4,5,5'-Hexachlor odiphenyl ether (13C12, 99%)	100	ng/mL