

## **PCN Native PAR Solution (unlabelled) 1000 ng/mL in Isooctane**

Art. ID	CIL-ECN-5497
Unit	1.2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
PCN 27	1,2,3,4-Tetrachloronaphthalene	[20020-02-4]	1000	ng/mL		
PCN 52	1,2,3,5,7-Pentachloronaphthalene	[53555-65-0]	1000	ng/mL		
PCN 66	1,2,3,4,6,7-Hexachloronaphthalene	[103426-96-6]	1000	ng/mL		
PCN 67	1,2,3,5,6,7-Hexachloronaphthalene	[103426-97-7]	1000	ng/mL		
PCN 68	1,2,3,5,6,8-Hexachloronaphthalene	[103426-95-5]	1000	ng/mL		
PCN 73	1,2,3,4,5,6,7-Heptachloronaphthalene	[58863-14-2]	1000	ng/mL		
PCN 75	Octachloronaphthalene	[2234-13-1]	100	ng/mL		