

EN-1948 Syringe Standard Solution

Art. ID CIL-ED-4140
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,3,4-Tetrachlorodibenzo-p-dioxin (13C12,99%)	[114423-99-3]	400	ng/mL		
	1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin (13C12,99%)	[109719-82-6]	400	ng/mL		