

## **Method 1613 Internal Standard Spiking Solution**

Art. ID CIL-EDF-5999  
Unit 0,5 mL  
Deliverydetails De Minimis UN1920 III E1

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,3,4-Tetrachlorodibenzo-p-dioxin (13C12,99%)	[114423-99-3]	200	ng/mL		
	1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin (13C12,99%)	[109719-82-6]	200	ng/mL		