

## **2,3,4,7-Tetrachlorodibenzofuran 50 µg/mL in Nonane**

|                 |   |
|-----------------|---|
| Art. ID         | CIL-EF-4031   |
| Unit            | 1,2 mL  |
| Deliverydetails | ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1 |

### Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter               | CAS number | Concentration/Value | Unit  | Method | Source |
|------------------|---------------------------------|------------|---------------------|-------|--------|--------|
|                  | 2,3,4,7-Tetrachlorodibenzofuran | [n/a]      | 50                  | µg/mL |        |        |