

Polybrominated Diphenyl Ether Performance Standard Solution

Art. ID CIL-EO-5101
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
BDE 77	3,3',4,4'-Tetrabromodiphenyl ether (13C12,99%) (BDE 77)	[93703-48-1]/un labelled	100	ng/mL		
BDE 126	3,3',4,4',5-Pentabromodiphenyl ether (13C12,99%) (BDE 126)	[1401355-56-3]	150	ng/mL		