

POPS Toxaphene surrogate solution with PCB syringe (13C10,99%)

Art. ID CIL-ES-5352-L
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------------------|------------|---------------------|-------|--------|--------|
| | trans-Chlordane (gamma) -13C10 | | 1 | µg/mL | | |