

PAH Native Standard Mixture (unlabelled)

Art. ID	CIL-ES-5438
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Naphthalene	[91-20-3]	200	µg/mL		
	Benz[a]anthracene	[56-55-3]	200	µg/mL		
	Phenanthrene	[85-01-8]	200	µg/mL		
	Fluoranthene	[206-44-0]	200	µg/mL		
	Benzo[b]fluoranthene	[205-99-2]	200	µg/mL		
	Benzo[a]pyrene	[50-32-8]	200	µg/mL		
	Benzo[g,h,i]perylene	[191-24-2]	200	µg/mL		
	Indeno[1,2,3-cd]pyrene	[193-39-5]	200	µg/mL		
	Dibenz[a,h]anthracene	[53-70-3]	200	µg/mL		
	Acenaphthylene	[208-96-8]	200	µg/mL		
	Acenaphthene	[83-32-9]	200	µg/mL		
	Fluorene	[86-73-7]	200	µg/mL		
	Pyrene	[129-00-0]	200	µg/mL		
	Benzo[k]fluoranthene	[207-08-9]	200	µg/mL		
	Perylene	[198-55-0]	200	µg/mL		
	Chrysene	[218-01-9]	200	µg/mL		