

### **Custom Solution PAH Mixture (D,98%) 2500 µg/mL in Toluene**

Art. ID CIL-ES-5481-5X5  
Unit 5 x 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

#### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chrysene (D12,98%)	[1719-03-5]	2500	µg/mL		
	Dibenz[a,h]anthracene (D14,98%)	[13250-98-1]	2500	µg/mL		
	Naphthalene (D8,99%)	[1146-65-2]	2500	µg/mL		
	Perylene (D12,98%)	[1520-96-3]	2500	µg/mL		
	Phenanthrene (D10,98%)	[1517-22-2]	2500	µg/mL		
	Acenaphthene (D10,98%)	[15067-26-2]	2500	µg/mL		