

PAH Two Component Mixture (D,98%) 2000 µg/mL in 80% isooctane/20% toluene

Art. ID CIL-ES-5498
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzo[a]pyrene (D12,98%)	[63466-71-7]	2000	ng/mL		
	Fluoranthene (D10,98%)	[93951-69-0]	2000	ng/mL		