

3,3',4,4'-Tetrabromobiphenyl (PBB-77) 100 µg/mL in Isooctane

Art. ID	CIL-PBB-77-CS
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,4'-Tetrabromobiphenyl (PBB-77)	[n/a]	100	µg/mL		