

## **2,3,3',4',6-Pentachlorobiphenyl 35 µg/mL in Isooctane**

|                 |                         |
|-----------------|-------------------------|
| Art. ID         | CIL-PCB-110             |
| Unit            | 1 mL                    |
| Deliverydetails | De Minimis UN1262 II E2 |

### Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter                             | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|---|--------------|---------------------|-------|--------|--------|
|                  | 2,3,3',4',6-Pentachloro<br>biphenyl (PCB 110) | [38380-03-9] | 35                  | µg/mL |        |        |